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The product range of DETARTRINE-Tube consists in ready-to-use prophylactic pastes for prophylaxis treatment after scaling and root planning, for the removal of drink, food and tobacco staining and for the polishing of dental surfaces.

DETARTRINE Product range

The abrasivity of a substance on the dentin is determined by the Relative Dentin Abrasion (RDA) index. The DETARTRINE product range comprises

3 products with different abrasivity to encompass every clinical situation:

- DETARTRINE 150Z (pink) - RDA 150 – Highly abrasive to efficiently remove heavy extrinsic stains.
- DETARTRINE 100ZF (white) – with Fluoride, RDA 100 – Moderately abrasive for most cleaning and polishing treatments, while participating in the remineralization of dental surfaces.
- DETARTRINE 50F (green) – with Fluoride, RDA 50 – Slightly abrasive, ideal for a gentle polishing of enamel and metallic, composite restorations, while participating in the remineralization of dental surfaces.

A full prophylaxis treatment begins with an abrasive paste and ends with a less abrasive one that contains fluoride.

With a purely mechanical abrasive action due to the presence of Zircon, DETARTRINE 150Z and DETARTRINE 100ZF pastes allow the removal of drink, food and tobacco stains and the polishing of the surface irregularities resulting from the use of manual or ultrasonic instruments, without eroding the enamel.

The silica contained in DETARTRINE 50F gives the paste polishing properties that allow the obtaining of smooth surfaces where dental plaque cannot easily adhere and a high lustre of the enamel, the metallic and composite restorations…

The associated fluoride contained in DETARTRINE 100ZF and 50F allows a superficial fluoridation of the enamel newly exposed by the scaling and root planning, protects from the chemical attacks and diminishes the hypersensitivity that often occurs after scaling.

PROFESSIONAL PROPHYLACTIC CLEANING PROCEDURE

The presence of plaque and tartar induces an inflammation with reversible or irreversible modifications of the periodontal tissues.

The regular and efficient removal of plaque and tartar by scaling is obtained with manual or ultrasonic instruments.

A final polishing is then necessary to delay the apposition of dental plaque and to remove the extrinsic drink, food and tobacco stains. This polishing is obtained with a brush covered with polishing paste, on a handpiece.

Scaling and polishing allow a thorough prophylactic cleaning of the tooth surfaces.

It is useful to complete this prophylaxis treatment with a conventional fluoridation treatment (gel, varnish, trays…).

INDICATIONS

Stain removal

Polishing after scaling and root planning

Polishing of restorations

CONTRAINDICATIONS

DETARTRINE prophylactic pastes are contraindicated on patients with a known allergy to one of the components of the paste.

INSTRUCTIONS FOR USE

Use a cup or a brush on a handpiece to apply the product.

Rinsing can be realized during or after the session.

Tightly close after use.

PRECAUTIONS FOR USE

Wear suitable protecting clothing and gloves.

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