# Septalloy

## NON GAMMA 2

NG 50 y NG 70

### Propiedades

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NG 50 and NG 70

### Properties

These new high copper dispersion amalgam systems are made of specially formulated micrograin powder of rounded particles. Prepared with mercury, the NON-GAMMA 2, NG 50 and NG 70 alloys make it possible to obtain an homogeneous, quickly setting material.

As most alloys of this type, SEPTALLOY NON-GAMMA 2, NG 50 and NG 70, show an increased resistance to corrosion as well as an exceptional dimensional stability.

### NG 50 - SEPTALLOY NON-GAMMA 2

#### Composition

Silver 50 %, Tin 30 % and Copper 20 %.

#### Directions for use

Alloy / mercury proportion  
- SEPTALLOY NON-GAMMA 2 - NG 50 ..... 1  
- MERCURY ..... 1

#### Amalgamation

With DENTOMAT or DUOMAT, allow about 30 seconds for mixing.

### NG 70 - SEPTALLOY NON-GAMMA 2

#### Composition

Ag. 70 %, Sn. 18.5 %, Cu. 11 % and Zn. 0.5 %.

#### Note

This alloy contains zinc and the amalgam made therefrom may exhibit corrosion and excessive expansion if contamination by moisture occurs during mixing or condensing. Consequently, the amalgam should neither be mixed with bare hands nor condensed into a wet cavity. Contamination by moisture must be avoided.

#### Directions for use

Alloy / mercury proportion  
- SEPTALLOY NON-GAMMA 2 - NG 70 ..... 1  
- MERCURY ..... 1,2

#### Amalgamation

With DENTOMAT or DUOMAT, allow about 30 seconds for mixing.

Physical characteristics	NG 50	NG 70
SEPTALLOY- MERCURY ratio	1= 1	1= 1,2
Mixing time in seconds	30	30
Working time in minutes	6 to 9	7 to 9
Compression strength in N/mm <sup>2</sup> after 1 hour after 24 hours	> 180 > 450	130 to 160 430 to 490
Flow in % (crushing under 10 N/mm <sup>2</sup> at 37°C for 3 to 24 hours after condensation)	< 0,3	< 0,3
Static Creep in % (crushing under 35 N/mm <sup>2</sup> at 37°C in 3 hours, 7 days after condensation)	< 0,2	< 0,2
Expansion in µm/m after 24 hours at 37°C	+ 3	+ 3

#### Presentations

##### NG 50 - SEPTALLOY NON-GAMMA 2

250 g bottle.

##### NG 70 - SEPTALLOY NON-GAMMA 2

250 g bottle.

**For professional use only in the practice of dentistry.**

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