COMPOSITION

DETARTRINE 150 Z :

Silicate de zirconium, silice, gomme adragante, arôme menthe, colorant, excipients.

DETARTRINE 100ZF :

Oxyde de zirconium, fluorure de sodium, gomme adragante, arôme menthe, excipients.

PROPRIÉTÉS

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COMPOSITION

DETARTRINE 150 Z :

Zirconium silicate, silica, gum tragacanth, spearmint flavour, colouring agent, excipients.

DETARTRINE 100ZF :

Zirconium oxide, sodium fluoride, gum tragacanth, spearmint flavour, excipients.

PROPERTIES

The product range of DETARTRINE-Tube consists in ready-to-use prophylactic pastes for prophylaxis treatment after scaling and root planning, for the removal of drink, food and tobacco staining and for the polishing of dental surfaces.

The abrasivity of a substance on the dentin is determined by the Relative Dentin Abrasion (RDA) index. The DETARTRINE product range comprises 2 products with different abrasivity to encompass every clinical situation:

- DETARTRINE 150Z (pink) - RDA 150 - Highly abrasive to efficiently remove heavy extrinsic stains.
- DETARTRINE 100ZF (white) - with Fluoride, RDA 100 - Moderately abrasive for most cleaning and polishing treatments, while participating in the remineralization of dental surfaces.

A full prophylaxis treatment begins with an abrasive paste and ends with a less abrasive one that contains fluoride.

With a purely mechanical abrasive action due to the presence of Zircon, DETARTRINE 150Z and DETARTRINE 100ZF pastes allow the removal of drink, food and tobacco stains and the polishing of the surface irregularities resulting from the use of manual or ultrasonic instruments, without eroding the enamel.

The associated fluoride contained in DETARTRINE 100ZF allows a superficial fluoridation of the enamel newly exposed by the scaling and root planning, protects from the chemical attacks and diminishes the hypersensitivity that often occurs after scaling.

INDICATIONS

Stain removal.
Polishing after scaling and root planning.
Polishing of restorations.

CONTRAINDICTION

Hypersensitivity to any component.

INSTRUCTIONS FOR USE

Professional Prophylactic Cleaning procedure

The presence of plaque and tartar induces an inflammation with reversible or irreversible modifications of the periodontal tissues.

The regular and efficient removal of plaque and tartar by scaling is obtained with manual or ultrasonic instruments.

A final polishing is then necessary to delay the apposition of dental plaque and to remove the extrinsic drink, food and tobacco stains. This polishing is obtained with a brush covered with polishing paste, on a handpiece.

Scaling and polishing allow a thorough prophylactic cleaning of the tooth surfaces.

It is useful to complete this prophylaxis treatment with a conventional fluoridation treatment (gel, varnish, trays...).

Instructions for use

- Use a cup or a brush on a handpiece to apply the product.
- Rinsing can be realized during or after the session.
- Tightly close after use.

WARNINGS AND PRECAUTIONS FOR USE

Wear suitable protecting clothing and gloves.

PRESENTATIONS

DETARTRINE 150Z :tube containing 65 mL of paste

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