DESCRIPTION

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DESCRIPTION

R.T.R. cone is made from β tricalcium phosphate granules of synthetic origin coated with a matrix of highly purified collagen fibres of bovine origin. R.T.R. calcium phosphates are manufactured using the latest technologies, which guarantees an extremely high purity level of the products.

PROPERTIES

R.T.R. is available as granules with a size ranging from 500 µm to 1 mm. The size of macropores varies from 100 µm to 400 µm and the micropores are less than 10 µm in diameter. These specific structural properties make it possible for newly formed bone to colonise macropores.

R.T.R. is gradually resorbable.

INDICATIONS

R.T.R. cone is designed for the filling and the reconstruction of bone defects in maxillofacial and dental surgery, after being impregnated with the patient's blood or with physiological solution (particularly for the filling of tooth sockets for instance).

R.T.R. cone is indicated particularly when the socket cannot be closed, since granules would probably leak out of the socket.

CONTRAINDICATIONS :

- Acute osteomyelitis or necrosis of the recipient site.
- Bone degeneration.
- Known active infection.
- Known allergy to any constituent of the bone substitute.

WARNINGS AND PRECAUTIONS FOR USE

The manufacturer of R.T.R. guarantees the material and manufacture quality.

Several criteria are of crucial importance for the operation with R.T.R. to be successful:

- appropriate selection of the patient who will undergo the operation,
- rigorous asepsis.

The practitioner is responsible for any complication which may result from improper use, faulty operative technique or lack of asepsis. These complications cannot be attributed to the manufacturer.

Underlying oral pathologies such as infections should be treated and ruled out before applying R.T.R. Avoid saliva contamination.

The bone site which will be in contact with R.T.R. should be debrided and, if possible, well vascularised, so as to create a high quality recipient bed.

Any non-used part of R.T.R. should be systematically disposed of since the product is no longer sterile.

Check the packaging integrity before use.

Do not use if packaging has been damaged or opened before first use.

Do not resterilise.

RESTERILISATION AND REUSE

The resterilisation and/or reuse of the R.T.R. bone substitutes is strictly prohibited, as it could expose the patient concerned to contamination and a subsequent risk of infection. Reusing a R.T.R. bone substitute would reduce its efficiency (decrease in mechanical strength).

INSTRUCTIONS FOR USE

R.T.R. cone is inserted into the dental socket after waiting that the socket has been filled with blood. It is always better to suture the operative site. However, should suturing not be possible, make sure that the cone has been completely covered with and protected by the wound edges.

In any case, the operative site will have to be closed by joining together the wound edges (coaptation) with suture stitches.

EXPECTED PERFORMANCE

As their chemical composition is very close to that of natural mineral bone, calcium phosphates – and particularly β-tricalcium phosphate – play an important role in the biological processes governing the stability and regeneration of bone tissue. Accordingly, the role of R.T.R. bone substitutes is to temporarily replace bone tissue in the event of major tissue deficiency caused by trauma or a physiological problem. Thus these devices enable initiation of the bone remodelling phenomenon whilst facilitating their replacement with newly formed natural bone.

STORAGE

Keep only in the original container.

EXPIRY DATE

Do not use after the expiry date mentioned on the outer package.

PRESENTATIONS

2 x 0,3 cm³ cone (Ø 6 mm, H 10 mm) in single-unit package sterilised by gamma radiation (minimal dose: 25 kGy).

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