Storage

Store between 5° and 27°C (41°-80°F)

Specifications

	ISO 4823 ADA 19	Oral setting time (min.)	Total setting time *	Strain in compression (min-max)	Shore-A hardness
Putty Hard	Type 0 Putty Consist.	3 min. 30 sec.	5 min. 30 sec.	1-3%	70+/-2
Putty Soft	Type 0 Putty Consist.	3 min. 30 sec.	5 min. 30 sec.	1-3%	60+/-2
Putty Soft Rapide	Type 0 Putty Consist.	2 min. 30 sec.	4 min.	1-3%	60+/-2
Putty Super Soft Rapide	Type 0 Putty Consist.	2 min. 30 sec.	4 min.	1-3%	60+/-2
Tray Material	Type 1 Heavy-bodied Consist.	2 min. 30 sec.	4 min.	3-5%	60+/-2
Regular	Type 2 Medium-bodied Consist.	3 min. 30 sec.	5 min. 30 sec.	3-5%	45+/-2
Light/Super Light	Type 3 Light-bodied Consist.	3 min. 30 sec.	5 min. 30 sec.	3-5%	45+/-2
Light Rapide	Type 3 Light-bodied Consist.	2 min. 30 sec.	4 min.	3-5%	45+/-2

* Time elapsed between the start of mixing and the final oral set, at a temperature of 23°C.

Dimensional changes < - 0.20 %

Changes in work time and temperature produce variations in setting time.

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- PERFEXIL PUTTY HARD : (grey) 1 x 280 ml jar of base
- PERFEXIL PUTTY SOFT : (green) 1 x 280 ml jar of catalyst
- PERFEXIL PUTTY SOFT RAPIDE : (blue)
- PERFEXIL PUTTY SUPERSOFT RAPIDE : (turquoise blue)
- PERFEXIL TRAY MATERIAL for dispensing gun 1 box containing 2 x 380 ml cartridges + 2 locking rings + 20 dynamic tips
- PERFEXIL CARTRIDGES TRAY MATERIAL : (khaki) for dispensing gun
- PERFEXIL CARTRIDGES LIGHT : (salmon-pink) 1 box containing 2 x 50 ml cartridges
- PERFEXIL CARTRIDGES REGULAR : (purple) + 12 mixing tips
- PERFEXIL CARTRIDGES LIGHT RAPIDE : (orange)
- PERFEXIL CARTRIDGES SUPER LIGHT : (yellow)
- DYNAMIC TIPS: bag of 50 pieces
- LOCKING RINGS : bag of 2 pieces
- MIXING TIPS: bag of 48 pieces
- INTRAORAL TIPS: bag of 96 pieces
- PERFEXIL ADHESIVE FOR RESIN TRAYS 5 ml bottle
- PERFEXIL UNIVERSAL ADHESIVE FOR STANDARD METAL TRAYS 10 ml bottle

For professional use only in the practice of dentistry.

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Perfexil platinum +

Addition silicone (vinyl polysiloxane)

New hydrophilic Vinyl Polysiloxane addition silicone.

The hydrophilic property of PERFEXIL makes it water-compatible and reduces the risk of bubbles.

The thixotropy of the material provides optimal flowability in the subgingival areas.

Indications

PERFEXIL is available in several viscosities, to accommodate the different impression techniques.

Instructions for use

- **Jars:** PUTTY HARD - PUTTY SOFT - PUTTY SOFT RAPIDE - PUTTY SUPERSOFT RAPIDE

Take equal amounts of both components (base and catalyst) using the relevant measuring spoons. Mix for about 30 seconds until you obtain a uniform, evenly-colored mixture with no streaks.

- **Cartridges for dispensing gun:** TRAY-MATERIAL - REGULAR - LIGHT - LIGHT RAPIDE - SUPER LIGHT

Assemble the automix device. Before attaching a tip, check that both components (base and catalyst) flow out of the cartridge evenly, by pressing gently on the dispenser lever. Extrude a small amount of product and discard it. Then attach the mixing tip to the cartridge and add an intraoral tip, if necessary. Close the jars immediately after use. Do not interchange caps and spoons.

- **Cartridges for dispensing machine (ratio 5 :1):** TRAY-MATERIAL

High viscosity, heavy-bodied consistency silicone, indicated for use with the simultaneous technique for loading of the tray, in combination with the PERFEXIL Platinum + range.

- After removing the cartridge from the package, pull the protection tab off the extrusion orifices and put it aside. You do not need to put too much force on the tab but hold the cartridge firmly so that you do not accidentally drop it, which may result in damage.

- Insert the cartridge into the dispensing machine (see manufacturer's instructions for use for the dispensing device).

- Before first use, you must ensure even extrusion of material from each side of the cartridge: start the machine until both components are evenly dispensed. Wipe off the extruded material.

- Insert the dynamic mixing tip. Then place the locking ring and turn it a quarter turn to lock the cartridge tip in place. The device is now ready for use.

- After use, it is recommended to remove the dynamic mixing tip and to place the plastic tab back on the cartridge. Discard the used mixing tip but keep the locking ring.

- To remove the cartridge, strictly follow the manufacturer's instructions for use of the machine. If the machine will not be used for several days, it is recommended to store the cartridge in an upright position.

- Important: always ensure even dispensing of the material before inserting the dynamic mixing tip. If you failed to do so, run the machine until an evenly-colored material is obtained, and discard the non-homogeneous mixture, before starting to use the product. Use exclusively the original locking rings and dynamic mixing tips supplied with TRAY MATERIAL. The use of other mixing tips will not ensure optimal mixing and may impair the quality of the product and adversely affect the final result. For further information about the use of the dispensing machine, please refer to the manufacturer's instructions.

Warnings and precautions for use

- Vinyl polysiloxanes may cause skin irritation or other allergic reactions in susceptible individuals.

- These silicones are chemically resistant, avoid contact with textiles.

Note: do not combine with condensation silicone.

Disinfection of impressions

The impressions made with PERFEXIL may be safely disinfected with suitable disinfecting solutions.

Casting

No need to wait, the impression may be poured immediately after disinfection and up to two weeks, provided that the impression is stored at room temperature.

The product is compatible with the best plasters on the market (type 4 and type 3) and with polyurethane resins.

Galvanization

The impressions can be silver or copper-plated safely.