fr

COMPOSITION

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en

COMPOSITION

Per sponge: native, non-denaturated, freeze-dried collagen of bovine origin..... 0.03g in average

TSE/BSE

The risk of virus transmission and animal transmissible spongiform encephalopathies has been assessed in compliance with the provisions of Regulation 722/2012.

INDICATIONS

Local haemostasis for dental use when control of bleeding by ligaturing or other conventional means is ineffective or impractical.

CONTRAINDICATION

Patients with a known history of allergy to collagen of bovine origin.

SIDE EFFECTS

Allergic reaction due to the presence of collagen of bovine origin.

INSTRUCTIONS FOR USE

Before inserting the haemostatic sponge, re-examine the socket. The inside of the socket should be thoroughly curetted to remove all granulation tissue and bone debris in order to reduce occurrence of alveolgia (dry socket). Apply a dry sponge with or without pressure, depending on the amount of bleeding, making sure that the socket does not dry out.

If necessary, the sponge can be held in place with a cross stitch over the socket.

The sponge will be completely and naturally resorbed.

WARNINGS AND PRECAUTIONS FOR USE

- Do not use in infected areas.
- Do not moisten the sponge.
- Any unused sponge remnants should be discarded, as no longer sterile.
- It is recommended to :
 - use dry, smooth instruments,
 - apply small pieces of sponge, just large enough to fill the socket.
- For single use only. Any reuse exposes other patients and the user to critical risks, including infection.
- Any excess not removed at the time of surgery may result in a foreign body reaction that may present with or without clinical signs and symptoms as a lesion or postoperative abscess formation upon imaging. Imaging may initially not be able to tell the difference.
- Check package integrity before use.

 **Do not use if package is open or damaged before first use.**

Do not re-sterilize.

STORAGE

To maintain sterility of the sponges and prevent moisture uptake, do not remove the heat-sealed film until use. Do not store above 25°C.

PRESENTATION

Box containing 3 blisters of 8 single use sponges. Sponges radio-sterilized by β-ray irradiation.

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