

Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 09.04.2018

Version number 2

Revision: 09.04.2018

SECTION 1: Identification of the substance/mixture and of the company/undertaking**· 1.1 Product identifier****· Trade name: Quit Janet enzym****· Article number / formula: N131****· 1.2 Relevant identified uses of the substance or mixture and uses advised against**

No further relevant information available.

· Application of the substance / the preparation : dental product**· 1.3 Details of the supplier of the safety data sheet****· Manufacturer/Supplier:**

SEPTODONT.

58 rue du Pont de Creteil

94100 Saint Maur des Fosses

France

Téléphone: +33(0)149767000

Fax: +33(0)148855401

msds@septodont.com

http://www.septodont.com

Made by :

Laboratoires ANIOS.

Adresse : PAVE DU MOULIN .59260.LILLE - HELLEMMES.FRANCE.

Téléphone : + 33 (0)3 20 67 67 67. Fax : + 33 (0)3 20 67 67 68.

e:mail : fds@anios.com

www.anios.com

· 1.4 Emergency telephone number:

France (INRS 24h/24) - Tel : +33 (0)1 45 42 59 59

Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112

Belgique - Tel : 32 070/245 245

UK : Call NHS 999 in a medical emergency. If you urgently need medical help or advice but it's not a life-threatening situation, call 111 (24/7). The NHS 111 number is currently only available in certain areas. If you're outside of these areas, you should call NHS Direct on 0845 4647

SECTION 2: Hazards identification**· 2.1 Classification of the substance or mixture****· Classification according to Regulation (EC) No 1272/2008 :**

GHS05

Skin Corr. 1B H314 Causes severe skin burns and eye damage.



GHS09

Aquatic Chronic 1 H410 Very toxic to aquatic life with long lasting effects.

· 2.2 Label elements**· Labelling according to Regulation (EC) No 1272/2008**

The product is classified and labelled according to the CLP regulation.

(Contd. on page 2)

GB

Safety data sheet

According to 1907/2006/EC, Article 31

Printing date 09.04.2018

Version number 2

Revision: 09.04.2018

Trade name: Quitanet enzym

(Contd. of page 1)

Hazard pictograms


GHS05 GHS09

Signal word Danger
Hazard-determining components of labelling:

 Didecylmethylpoly(oxethyl) Ammonium Propionate
 AMINES,N-C12-14-ALKYLTRIMETHYLENEDI-

Hazard statements

 H314 Causes severe skin burns and eye damage.
 H410 Very toxic to aquatic life with long lasting effects.

Precautionary statements

P102 Keep out of reach of children.
 P273 Avoid release to the environment.
 P280 Wear protective gloves/protective clothing/eye protection/face protection.
 P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
 P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower].
 P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 P310 Immediately call a POISON CENTER/doctor.
 P501 Dispose of contents/container in accordance with local/regional/national/international regulations.





2.3 Other hazards
Results of PBT and vPvB assessment

- **PBT:** Not applicable.
- **vPvB:** Not applicable.

SECTION 3: Composition/information on ingredients

3.2 Chemical characterisation: Mixtures
Description: Mixture: consisting of the following components.

Dangerous components:

CAS: 94667-33-1	Didecylmethylpoly(oxethyl) Ammonium Propionate 	10-25%
CAS: 71060-57-6	Alkyl polyglycoléther C8-10 (OE) 	5-10%
CAS: 90640-43-0 Reg.nr.: 01-2119957843-25	AMINES,N-C12-14-ALKYLTRIMETHYLENEDI- 	≤2.5%
CAS: 18472-51-0 EINECS: 242-354-0 Reg.nr.: 01-2119946568-22-0001	D-gluconic acid, compound with N,N''-bis(4-chlorophenyl)-3,12-diimino-2,4,11,13-tetraazatetradecanediamidine (2:1) 	≤1%

SVHC No

(Contd. on page 3)

Safety data sheet

According to 1907/2006/EC, Article 31

Printing date 09.04.2018

Version number 2

Revision: 09.04.2018

Trade name: Quitanet enzym

(Contd. of page 2)

- **Additional information:** For the wording of the listed risk phrases, refer to section 16.

SECTION 4: First aid measures

· 4.1 Description of first aid measures

· General information:

Immediately remove any clothing soiled by the product.

Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.

· After excessive inhalation:

In case of unconsciousness place patient in a stable laying down side position for transportation.

· After skin contact: Immediately wash with water and soap and rinse thoroughly.

· After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.

· After swallowing:

Do not induce vomiting; call for medical help immediately.

Call for a doctor immediately.

Drink plenty of water and provide fresh air. Call for a doctor immediately.

· 4.2 Most important symptoms and effects, both acute and delayed :

No further relevant information available.

· 4.3 Indication of any immediate medical attention and special treatment needed :

No further relevant information available.

SECTION 5: Firefighting measures

· 5.1 Suitable extinguishing agents:

CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

· 5.2 Special hazards arising from the substance or mixture : No further relevant information available.

· 5.3 Advice for firefighters

· Protective equipment: No special measures required.

SECTION 6: Accidental release measures

· 6.1 Personal precautions, protective equipment and emergency procedures :

Wear protective equipment. Keep unprotected persons away.

· 6.2 Environmental precautions:

Do not allow product to reach sewage system or any water course.

Inform respective authorities in case of seepage into water course or sewage system.

Dilute with plenty of water.

Do not allow the product to enter sewers/ surface water or ground water.

· 6.3 Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Use neutralising agent.

Dispose contaminated material as waste according to item 13.

Ensure adequate ventilation.

· 6.4 Reference to other sections :

See Section 7 for information on safe handling.

See Section 8 for information on personal protective equipment.

See Section 13 for disposal information.

GB

(Contd. on page 4)

Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 09.04.2018

Version number 2

Revision: 09.04.2018

Trade name: Quitanet enzym

(Contd. of page 3)

SECTION 7: Handling and storage

- **7.1 Precautions for safe handling :**
Ensure good ventilation/exhaustion at the workplace.
Prevent formation of aerosols.
- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**
Protect from frost.
Keep container tightly sealed.
- **Recommended storage temperature:** 5°C-35°C
- **7.3 Specific end use(s) :** No further relevant information available.

SECTION 8: Exposure controls/personal protection

- **Additional information about design of technical facilities:** No further data; see section 7.

· **8.1 Control parameters**

- **Ingredients with limit values that require monitoring at the workplace:**

CAS: 107-21-1 ethanediol (<1%)

WEL (Great Britain)	Short-term value: 104** mg/m ³ , 40** ppm Long-term value: 10* 52** mg/m ³ , 20** ppm Sk *particulate **vapour
IOELV (EU)	Short-term value: 104 mg/m ³ , 40 ppm Long-term value: 52 mg/m ³ , 20 ppm Skin
MAK (Switzerland)	Short-term value: 52 mg/m ³ , 20 ppm Long-term value: 26 mg/m ³ , 10 ppm H SSc;

- **Additional information:**
The most current valid lists have been used as a basis for the production of this document.

· **8.2 Exposure controls**

- **Personal protective equipment:**
- **General protective and hygienic measures:**
Keep away from food, beverages and petfood.
Immediately remove all soiled and contaminated clothing
Wash hands before breaks and at the end of work.
Avoid contact with the eyes and skin.
- **Respiratory protection:** Not required.
- **Protection of hands:**



Protective gloves

- **Gloves material :** Butyl rubber, BR

(Contd. on page 5)

Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 09.04.2018

Version number 2

Revision: 09.04.2018

Trade name: Quitanet enzym

· **Eye protection:**

(Contd. of page 4)



Tightly sealed goggles

SECTION 9: Physical and chemical properties

· **9.1 Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· Form:	Liquid
· Colour:	Blue
· Odour:	Characteristic
· Odour threshold:	Not determined.

· **pH-value:** 6-9.5

· **Change in condition**

· Melting point/freezing point:	Undetermined.
· Initial boiling point and boiling range:	Undetermined.

· **Flash point:** Not applicable.

· **Flammability (solid, gas):** Not applicable.

· **Ignition temperature:** 425 °C

· **Decomposition temperature:** Not determined.

· **Auto-ignition temperature:** Product is not self igniting.

· **Explosive properties:** Product does not present an explosion hazard.

· **Explosion limits:**

· Lower:	Not determined.
· Upper:	Not determined.

· **Vapour pressure:** Not determined.

· Density at 20 °C:	1 g/cm ³
· Relative density :	Not determined.
· Vapour density :	Not determined.
· Evaporation rate :	Not determined.

· **Solubility in / Miscibility with water:** Fully miscible.

· **Segregation coefficient (n-octanol/water):** Not determined.

· **Viscosity:**

· Dynamic:	Not determined.
· Kinematic:	Not determined.

· **9.2 Other information :** No further relevant information available.

Safety data sheet

According to 1907/2006/EC, Article 31

Printing date 09.04.2018

Version number 2

Revision: 09.04.2018

 Trade name: *Quit Janet enzym*

(Contd. of page 5)

SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- **10.4 Conditions to avoid :** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

SECTION 11: Toxicological information

- **11.1 Information on toxicological effects**
- **Acute toxicity** Based on available data, the classification criteria are not met.

- **LD/LC50 values relevant for classification:**

CAS: 94667-33-1 Didecylmethylpoly(oxethyl) Ammonium Propionate

Oral	LD50	500 mg/kg (ATE)
------	------	-----------------

CAS: 71060-57-6 Alkyl polyglycoléther C8-10 (OE)

Oral	LD50	500 mg/kg (ATE)
------	------	-----------------

CAS: 90640-43-0 AMINES,N-C12-14-ALKYLTRIMETHYLENEDI-

Oral	LD50	100 mg/kg (ATE)
------	------	-----------------

CAS: 18472-51-0 D-gluconic acid, compound with N,N''-bis(4-chlorophenyl)-3,12-diimino-2,4,11,13-tetraazatetradecanediamidine (2:1)

Oral	LD50	>2,000 mg/kg (rat)
------	------	--------------------

Dermal	LD50	>2,000 mg/kg (rabbit)
--------	------	-----------------------

- **Primary irritant effect:**
- **On the skin:**
Causes severe skin burns and eye damage.
- **On the eyes:**
Causes severe skin burns and eye damage.
- **Respiratory or skin sensitization** Based on available data, the classification criteria are not met.
- **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

- **12.1 Toxicity**

- **Aquatic toxicity:**

CAS: 18472-51-0 D-gluconic acid, compound with N,N''-bis(4-chlorophenyl)-3,12-diimino-2,4,11,13-tetraazatetradecanediamidine (2:1)

LC50 96h	2.08 mg/l (fish)
----------	------------------

EC50 72h	0.081 mg/l (Algae)
----------	--------------------

(Contd. on page 7)

Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 09.04.2018

Version number 2

Revision: 09.04.2018

Trade name: Quitanet enzym

(Contd. of page 6)

EC50 48h 0.087 mg/l (Crustaceans)

- **12.2 Persistence and degradability** : No further relevant information available.
- **12.3 Bioaccumulative potential** : No further relevant information available.
- **12.4 Mobility in soil** : No further relevant information available.
- **Ecotoxicological effects:**
- **Remark:** Very toxic for fish
- **Additional ecological information:**
- **General notes:**
Water hazard class 3 (German Regulation) (Self-assessment): extremely hazardous for water
Do not allow product to reach ground water, water course or sewage system, even in small quantities.
Must not reach sewage water or drainage ditch undiluted or unneutralised.
Danger to drinking water if even extremely small quantities leak into the ground.
Also poisonous for fish and plankton in water bodies.
Very toxic for aquatic organisms
- **12.5 Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **12.6 Other adverse effects** : No further relevant information available.

SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation :**



Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

· **European waste catalogue**

18 00 00	WASTES FROM HUMAN OR ANIMAL HEALTH CARE AND/OR RELATED RESEARCH (EXCEPT KITCHEN AND RESTAURANT WASTES NOT ARISING FROM IMMEDIATE HEALTH CARE)
18 01 00	wastes from natal care, diagnosis, treatment or prevention of disease in humans
18 01 06*	chemicals consisting of or containing hazardous substances
HP 5	Specific Target Organ Toxicity (STOT)/Aspiration Toxicity
HP 6	Acute Toxicity
HP 8	Corrosive
HP 14	Ecotoxic

- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.
- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

SECTION 14: Transport information

- **14.1 UN-Number**
- **ADR, IMDG, IATA** UN1903

(Contd. on page 8)

Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 09.04.2018

Version number 2

Revision: 09.04.2018

Trade name: Quitanet enzym

(Contd. of page 7)

· 14.2 UN proper shipping name

· ADR

UN1903 DISINFECTANT, LIQUID, CORROSIVE, N.O.S.
(Didecylmethylpoly(oxethyl) Ammonium Propionate,
AMINES,N-C12-14-ALKYLTRIMETHYLENEDI-),
ENVIRONMENTALLY HAZARDOUS

· IMDG

DISINFECTANT, LIQUID, CORROSIVE, N.O.S.
(Didecylmethylpoly(oxethyl) Ammonium Propionate,
AMINES,N-C12-14-ALKYLTRIMETHYLENEDI-, D-
gluconic acid, compound with N,N''-bis(4-chlorophenyl)-
3,12-diimino-2,4,11,13-tetraazatetradecanediamidine (2:
1)), MARINE POLLUTANT

· IATA

DISINFECTANT, LIQUID, CORROSIVE, N.O.S.
(Didecylmethylpoly(oxethyl) Ammonium Propionate,
AMINES,N-C12-14-ALKYLTRIMETHYLENEDI-)

· 14.3 Transport hazard class(es)

· ADR



· Class

8 (C9) Corrosive substances.

· Label

8

· IMDG



· Class

8 Corrosive substances.

· Label

8

· IATA



· Class

8 Corrosive substances.

· Label

8

· 14.4 Packing group

· ADR, IMDG, IATA

II

· 14.5 Environmental hazards:

Product contains environmentally hazardous substances:
Didecylmethylpoly(oxethyl) Ammonium Propionate

· Marine pollutant:

Yes

· Special marking (ADR):

Symbol (fish and tree)
Symbol (fish and tree)

· 14.6 Special precautions for user

Warning: Corrosive substances.

· Danger code (Kemler):

86

· EMS Number:

F-A,S-B

(Contd. on page 9)

Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 09.04.2018

Version number 2

Revision: 09.04.2018

Trade name: Quitanet enzym

(Contd. of page 8)

· Stowage Category	B
· 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code	<i>Not applicable.</i>
· Transport/Additional information:	
· ADR	
· Limited quantities (LQ)	1L
· Excepted quantities (EQ)	Code: E2 Maximum net quantity per inner packaging: 30 ml Maximum net quantity per outer packaging: 500 ml
· Transport category	2
· Tunnel restriction code	E
· IMDG	
· Limited quantities (LQ)	1L
· Excepted quantities (EQ)	Code: E2 Maximum net quantity per inner packaging: 30 ml Maximum net quantity per outer packaging: 500 ml
· UN "Model Regulation":	UN 1903 DISINFECTANT, LIQUID, CORROSIVE, N.O.S. (DIDECYLMETHYLPOLY(OXETHYL) AMMONIUM PROPIONATE, AMINES, N-C12-14-ALKYLTRIMETHYLENEDI-), 8, II, ENVIRONMENTALLY HAZARDOUS

SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **Seveso category E1** Hazardous to the Aquatic Environment
- **Qualifying quantity (tonnes) for the application of lower-tier requirements** 100 t
- **Qualifying quantity (tonnes) for the application of upper-tier requirements** 200 t
- **REGULATION (EC) No 1907/2006 ANNEX XVII** Conditions of restriction: 3
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Relevant phrases**
- H301 Toxic if swallowed.
- H302 Harmful if swallowed.
- H314 Causes severe skin burns and eye damage.
- H318 Causes serious eye damage.
- H372 Causes damage to organs through prolonged or repeated exposure.
- H400 Very toxic to aquatic life.
- H410 Very toxic to aquatic life with long lasting effects.

(Contd. on page 10)

Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 09.04.2018

Version number 2

Revision: 09.04.2018

Trade name: Quitanet enzym

(Contd. of page 9)

· Classification according to Regulation (EC) No 1272/2008

Skin corrosion/irritation
Hazardous to the aquatic environment - long-term
(chronic) aquatic hazard

The classification of the mixture is generally based on the calculation method using substance data according to Regulation (EC) No 1272/2008.

· Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
IMDG: International Maritime Code for Dangerous Goods
IATA: International Air Transport Association
GHS: Globally Harmonised System of Classification and Labelling of Chemicals
EINECS: European Inventory of Existing Commercial Chemical Substances
ELINCS: European List of Notified Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
LC50: Lethal concentration, 50 percent
LD50: Lethal dose, 50 percent
PBT: Persistent, Bioaccumulative and Toxic
SVHC: Substances of Very High Concern
vPvB: very Persistent and very Bioaccumulative
Acute Tox. 3: Acute toxicity – Category 3
Acute Tox. 4: Acute toxicity – Category 4
Skin Corr. 1B: Skin corrosion/irritation – Category 1B
Eye Dam. 1: Serious eye damage/eye irritation – Category 1
STOT RE 1: Specific target organ toxicity (repeated exposure) – Category 1
Aquatic Acute 1: Hazardous to the aquatic environment - acute aquatic hazard – Category 1
Aquatic Chronic 1: Hazardous to the aquatic environment - long-term aquatic hazard – Category 1