



**Safety data sheet**  
**According to 1907/2006/EC, Article 31**

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**SECTION 1: Identification of the substance/mixture and of the company/undertaking**

- **1.1 Product identifier**
- **Trade name:** CALYPSO ORANGE
- **Article number / formula:** 265S
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the preparation :** dental product
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
SEPTODONT.  
58 rue du Pont de Creteil  
94100 Saint Maur des Fosses  
France  
Téléphone: +33(0)149767000  
Fax: +33(0)148855401  
msds@septodont.com  
http://www.septodont.com
- **1.4 Emergency telephone number:**  
France (INRS 24h/24) - Tel : +33 (0)1 45 42 59 59  
Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112  
Belgique - Tel : 32 070/245 245  
UK : +44 (0)-870-190-6777 (CARECHEM24)

**SECTION 2: Hazards identification**

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008 :**  
The product is not classified according to the CLP regulation.
- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008** Not classified as hazardous
- **Hazard pictograms** Not classified as hazardous
- **Signal word** Not classified as hazardous
- **Hazard statements** Not classified as hazardous
- **2.3 Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

**SECTION 3: Composition/information on ingredients**

- **3.2 Chemical characterisation: Mixtures**
- **Description:** Mixture: consisting of the following components.
- **Dangerous components:** Not classified as hazardous
- **SVHC** No
- **Additional information:** For the wording of the listed risk phrases, refer to section 16.

**SECTION 4: First aid measures**

- **4.1 Description of first aid measures**
- **General information:** No special measures required.

(Contd. on page 2)

GB

## Safety data sheet

### According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: CALYPSO ORANGE**

(Contd. of page 1)

- **After excessive inhalation:** Provide fresh air; consult a doctor in case of medical complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** If symptoms persist consult a doctor.
- **4.2 Most important symptoms and effects, both acute and delayed :**  
No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed :**  
No further relevant information available.

### SECTION 5: Firefighting measures

- **5.1 Suitable extinguishing agents:** Use fire extinguishing methods suitable to surrounding conditions.
- **5.2 Special hazards arising from the substance or mixture :** No further relevant information available.
- **5.3 Advice for firefighters**
- **Protective equipment:** No special measures required.

### SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures :** Not required.
- **6.2 Environmental precautions:**  
Dilute with plenty of water.  
Do not allow the product to enter sewers/ surface water or ground water.
- **6.3 Methods and material for containment and cleaning up:**  
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **6.4 Reference to other sections :**  
No dangerous substances are released.  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protective equipment.  
See Section 13 for disposal information.

### SECTION 7: Handling and storage

- **7.1 Precautions for safe handling :** No special measures required.
- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **7.3 Specific end use(s) :** No further relevant information available.

### SECTION 8: Exposure controls/personal protection

- **Additional information about design of technical facilities:** No further data; see section 7.

(Contd. on page 3)

**Safety data sheet**  
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: CALYPSO ORANGE**

(Contd. of page 2)

**· 8.1 Control parameters**

**· Ingredients with limit values that require monitoring at the workplace:**

**CAS: 56-81-5 glycerol (50-100%)**

WEL (Great Britain)	Long-term value: 10 mg/m <sup>3</sup>
MAK (Switzerland)	Short-term value: 100 e mg/m <sup>3</sup> Long-term value: 50 e mg/m <sup>3</sup> SSc;

**CAS: 57-55-6 Methyl glycol (25-50%)**

WEL (Great Britain)	Long-term value: 474* 10** mg/m <sup>3</sup> , 150* ppm *total vapour and particulates **particulates
---------------------	--

**CAS: 64-17-5 ethanol (<1%)**

WEL (Great Britain)	Long-term value: 1920 mg/m <sup>3</sup> , 1000 ppm
MAK (Switzerland)	Short-term value: 1920 mg/m <sup>3</sup> , 1000 ppm Long-term value: 960 mg/m <sup>3</sup> , 500 ppm SSc;

**· Additional information:**

The most current valid lists have been used as a basis for the production of this document.

**· 8.2 Exposure controls**

**· Personal protective equipment:**

**· General protective and hygienic measures:**

The usual precautionary measures are to be adhered to when handling chemicals.

**· Respiratory protection:** Not required.

**· Protection of hands:** Not required.

**· Gloves material :** -

**· Eye protection:** Goggles recommended during refilling

**SECTION 9: Physical and chemical properties**

**· 9.1 Information on basic physical and chemical properties**

**· General Information**

**· Appearance:**

<b>Form:</b>	Liquid
<b>Colour:</b>	According to product specification
<b>· Odour:</b>	Characteristic
<b>· Odour threshold:</b>	Not determined.

**· pH-value:** Not determined.

**· Change in condition**

<b>Melting point/freezing point:</b>	Undetermined.
<b>Initial boiling point and boiling range:</b>	Undetermined.

**· Flash point:** Not applicable.

**· Flammability (solid, gas):** Not applicable.

**· Ignition temperature:** Not applicable.

**· Decomposition temperature:** Not determined.

**· Auto-ignition temperature:** Product is not self igniting.

(Contd. on page 4)

**Safety data sheet**  
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: CALYPSO ORANGE**

(Contd. of page 3)

· <b>Explosive properties:</b>	<i>Product does not present an explosion hazard.</i>
· <b>Explosion limits:</b>	
<b>Lower:</b>	<i>Not determined.</i>
<b>Upper:</b>	<i>Not determined.</i>
· <b>Vapour pressure at 20 °C:</b>	<i>0,1 hPa</i>
· <b>Density:</b>	<i>Not determined.</i>
· <b>Relative density :</b>	<i>Not determined.</i>
· <b>Vapour density :</b>	<i>Not determined.</i>
· <b>Evaporation rate :</b>	<i>Not determined.</i>
· <b>Solubility in / Miscibility with water:</b>	<i>Fully miscible.</i>
· <b>Segregation coefficient (n-octanol/water):</b>	<i>Not determined.</i>
· <b>Viscosity:</b>	
<b>Dynamic:</b>	<i>Not determined.</i>
<b>Kinematic:</b>	<i>Not determined.</i>
· <b>Solvent content:</b>	
<b>VOC (EC)</b>	<i>38.61 %</i>
· <b>9.2 Other information :</b>	<i>No further relevant information available.</i>

**SECTION 10: Stability and reactivity**

- **10.1 Reactivity** *No further relevant information available.*
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:** *No decomposition if used according to specifications.*
- **10.3 Possibility of hazardous reactions :** *No dangerous reactions known.*
- **10.4 Conditions to avoid :** *No further relevant information available.*
- **10.5 Incompatible materials:** *No further relevant information available.*
- **10.6 Hazardous decomposition products:** *No dangerous decomposition products known.*

**SECTION 11: Toxicological information**

- **11.1 Information on toxicological effects**
- **Acute toxicity** *Based on available data, the classification criteria are not met.*
- **Primary irritant effect:**
- **On the skin:** *Based on available data, the classification criteria are not met.*
- **On the eyes:** *Based on available data, the classification criteria are not met.*
- **Respiratory or skin sensitization** *Based on available data, the classification criteria are not met.*
- **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**
- **Germ cell mutagenicity** *Based on available data, the classification criteria are not met.*
- **Carcinogenicity** *Based on available data, the classification criteria are not met.*
- **Reproductive toxicity** *Based on available data, the classification criteria are not met.*
- **STOT-single exposure** *Based on available data, the classification criteria are not met.*
- **STOT-repeated exposure** *Based on available data, the classification criteria are not met.*
- **Aspiration hazard** *Based on available data, the classification criteria are not met.*

GB

(Contd. on page 5)



**Safety data sheet**  
**According to 1907/2006/EC, Article 31**

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: CALYPSO ORANGE**

(Contd. of page 4)

**SECTION 12: Ecological information**

- **12.1 Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability :** No further relevant information available.
- **12.3 Bioaccumulative potential :** No further relevant information available.
- **12.4 Mobility in soil :** No further relevant information available.
- **Additional ecological information:**
- **General notes:**  
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water.  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **12.5 Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **12.6 Other adverse effects :** No further relevant information available.

**SECTION 13: Disposal considerations**

- **13.1 Waste treatment methods**
- **Recommendation :** Smaller quantities can be disposed of with household waste.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.
- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

**SECTION 14: Transport information**

- |   |   |
|---|---|
| · <b>14.1 UN-Number</b><br>· <b>ADR, ADN, IMDG, IATA</b>                                    | Not classified as hazardous for transport |
| · <b>14.2 UN proper shipping name</b><br>· <b>ADR, ADN, IMDG, IATA</b>                      | Not classified as hazardous for transport |
| · <b>14.3 Transport hazard class(es)</b><br>· <b>ADR, ADN, IMDG, IATA</b><br>· <b>Class</b> | Not classified as hazardous for transport |
| · <b>14.4 Packing group</b><br>· <b>ADR, IMDG, IATA</b>                                     | Not classified as hazardous for transport |
| · <b>14.5 Environmental hazards:</b><br>· <b>Marine pollutant:</b>                          | No  |
| · <b>14.6 Special precautions for user</b>  | Not applicable.                           |
| · <b>14.7 Transport in bulk according to Annex II of Marpol and the IBC Code</b>            | Not applicable.                           |
| · <b>UN "Model Regulation":</b>   | Not classified as hazardous for transport |

-GB-

(Contd. on page 6)



**Safety data sheet**  
**According to 1907/2006/EC, Article 31**

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

**Trade name: CALYPSO ORANGE**

(Contd. of page 5)

**SECTION 15: Regulatory information**

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

**SECTION 16: Other information**

*This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*

**Abbreviations and acronyms:**

*ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)*  
*IMDG: International Maritime Code for Dangerous Goods*  
*IATA: International Air Transport Association*  
*GHS: Globally Harmonised System of Classification and Labelling of Chemicals*  
*EINECS: European Inventory of Existing Commercial Chemical Substances*  
*ELINCS: European List of Notified Chemical Substances*  
*CAS: Chemical Abstracts Service (division of the American Chemical Society)*  
*VOC: Volatile Organic Compounds (USA, EU)*  
*PBT: Persistent, Bioaccumulative and Toxic*  
*SVHC: Substances of Very High Concern*  
*vPvB: very Persistent and very Bioaccumulative*