



Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

SECTION 1: Identification of the substance/mixture and of the company/undertaking

- **1.1 Product identifier**
- **Trade name: DETARTRINE 150 Z**
- **Article number / formula: 274S**
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
No further relevant information available.
- **Application of the substance / the preparation : dental product**
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
SEPTODONT.
58 rue du Pont de Creteil
94100 Saint Maur des Fosses
France
Téléphone: +33(0)149767000
Fax: +33(0)148855401
msds@septodont.com
http://www.septodont.com
- **1.4 Emergency telephone number:**
France (INRS 24h/24) - Tel : +33 (0)1 45 42 59 59
Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112
Belgique - Tel : 32 070/245 245
UK : +44 (0)-870-190-6777 (CARECHEM24)

SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008 :**
The product is not classified according to the CLP regulation.
- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008** Not classified as hazardous
- **Hazard pictograms** Not classified as hazardous
- **Signal word** Not classified as hazardous
- **Hazard statements** Not classified as hazardous
- **Additional information:**
EUH208 Contains Menthone. May produce an allergic reaction.
EUH210 Safety data sheet available on request.
- **2.3 Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

SECTION 3: Composition/information on ingredients

- **3.2 Chemical characterisation: Mixtures**
- **Description:** Mixture: consisting of the following components.

(Contd. on page 2)

GB

Safety data sheet

According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

Trade name: DETARTRINE 150 Z

(Contd. of page 1)

Dangerous components:

CAS: 10101-52-7	Zirconium(IV) silicate ⚠ Acute Tox. 4, H332; Skin Irrit. 2, H315; Eye Irrit. 2, H319; STOT SE 3, H335	10-25%
CAS: 64-17-5 EINECS: 200-578-6 Index number: 603-002-00-5 RTECS: KQ 6300000 Reg.nr.: 01-2119457610-43-xxxx	ethanol ⚠ Flam. Liq. 2, H225; ⚠ Eye Irrit. 2, H319	≤2.5%

· SVHC No

 · **Additional information:** For the wording of the listed risk phrases, refer to section 16.

SECTION 4: First aid measures

- **4.1 Description of first aid measures**
- **General information:** No special measures required.
- **After excessive inhalation:** Provide fresh air; consult a doctor in case of medical complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** If symptoms persist consult a doctor.
- **4.2 Most important symptoms and effects, both acute and delayed :**
No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed :**
No further relevant information available.

SECTION 5: Firefighting measures

- **5.1 Suitable extinguishing agents:**
CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **5.2 Special hazards arising from the substance or mixture :** No further relevant information available.
- **5.3 Advice for firefighters**
- **Protective equipment:** No special measures required.

SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures :** Not required.
- **6.2 Environmental precautions:** Do not allow the product to enter sewers/ surface water or ground water.
- **6.3 Methods and material for containment and cleaning up:** Pick up mechanically.
- **6.4 Reference to other sections :**
No dangerous substances are released.
See Section 7 for information on safe handling.
See Section 8 for information on personal protective equipment.
See Section 13 for disposal information.

SECTION 7: Handling and storage

- **7.1 Precautions for safe handling :** No special measures required.
- **Information about fire - and explosion protection:** No special measures required.

(Contd. on page 3)

Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

Trade name: DETARTRINE 150 Z

(Contd. of page 2)

- 7.2 Conditions for safe storage, including any incompatibilities
- Storage:
- Requirements to be met by storerooms and receptacles: No special requirements.
- Information about storage in one common storage facility: Not required.
- Further information about storage conditions: None.
- 7.3 Specific end use(s) : No further relevant information available.

SECTION 8: Exposure controls/personal protection

- Additional information about design of technical facilities: No further data; see section 7.

· 8.1 Control parameters

- Ingredients with limit values that require monitoring at the workplace:

CAS: 14808-60-7 Quartz (SiO₂) (25-50%)

MAK (Switzerland)	Long-term value: 0.15 a mg/m ³ P C1a SSc;
-------------------	---

CAS: 56-81-5 glycerol (10-25%)

WEL (Great Britain)	Long-term value: 10 mg/m ³
MAK (Switzerland)	Short-term value: 100 e mg/m ³ Long-term value: 50 e mg/m ³ SSc;

CAS: 64-17-5 ethanol (≤2.5%)

WEL (Great Britain)	Long-term value: 1920 mg/m ³ , 1000 ppm
MAK (Switzerland)	Short-term value: 1920 mg/m ³ , 1000 ppm Long-term value: 960 mg/m ³ , 500 ppm SSc;

CAS: 1310-58-3 potassium hydroxide (<0.1%)

WEL (Great Britain)	Short-term value: 2 mg/m ³
MAK (Switzerland)	Long-term value: 2 e mg/m ³

- DNELs

CAS: 64-17-5 ethanol

Dermal	DNEL Worker long-term systemic	343 mg/kg bw/d
Inhalative	DNEL Worker long-term systemic	950 mg/m ³
	DNEL Worker acute/short-term local effects	19 mg/m ³

- PNECs

CAS: 64-17-5 ethanol

Oral	PNEC sediment (fresh water)	3.6 mg/kg
	PNEC STP	0.72 mg/L
	PNEC aqua (fresh water)	0.96 mg/L
	PNEC aqua (marine water)	0.79 mg/L
	PNEC sediment (marine water)	2.9 mg/kg
	PNEC soil	0.63 mg/kg

- Additional information:

The most current valid lists have been used as a basis for the production of this document.

(Contd. on page 4)

Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

Trade name: DETARTRINE 150 Z

(Contd. of page 3)

- **8.2 Exposure controls**
- **Personal protective equipment:**
- **General protective and hygienic measures:**
The usual precautionary measures are to be adhered to when handling chemicals.
- **Respiratory protection:** Not required.
- **Protection of hands:** Not required.
- **Gloves material :** -
- **Eye protection:** Not required.

SECTION 9: Physical and chemical properties

· **9.1 Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

- | | |
|---------------------------|-----------------|
| · Form: | Pasty |
| · Colour: | Pink |
| · Odour: | Characteristic |
| · Odour threshold: | Not determined. |

· **pH-value:** Not applicable.

· **Change in condition**

- | | |
|---|---------------|
| · Melting point/freezing point: | Undetermined. |
| · Initial boiling point and boiling range: | Undetermined. |

· **Flash point:** Not applicable.

· **Flammability (solid, gas):** Not determined.

· **Ignition temperature:** 400 °C

· **Decomposition temperature:** Not determined.

· **Auto-ignition temperature:** Product is not self igniting.

· **Explosive properties:** Product does not present an explosion hazard.

· **Explosion limits:**

- | | |
|-----------------|-----------------|
| · Lower: | 0,9 Vol % |
| · Upper: | Not determined. |

· **Vapour pressure at 20 °C:** 0,1 hPa

· **Density:** Not determined.

· **Relative density :** Not determined.

· **Vapour density :** Not applicable.

· **Evaporation rate :** Not applicable.

· **Solubility in / Miscibility with water:**

Insoluble.

· **Segregation coefficient (n-octanol/water):** Not determined.

· **Viscosity:**

- | | |
|---------------------|-----------------|
| · Dynamic: | Not applicable. |
| · Kinematic: | Not applicable. |

· **Solvent content:**

· **VOC (EC)** 0.00 %

(Contd. on page 5)

Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

Trade name: DETARTRINE 150 Z

(Contd. of page 4)

· **9.2 Other information :** No further relevant information available.

SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions :** No dangerous reactions known.
- **10.4 Conditions to avoid :** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

SECTION 11: Toxicological information

- **11.1 Information on toxicological effects**
- **Acute toxicity** Based on available data, the classification criteria are not met.

· **LD/LC50 values relevant for classification:**

CAS: 10101-52-7 Zirconium(IV) silicate

Inhalative	LC50/4 h	1.5 mg/l (ATE)
------------	----------	----------------

CAS: 64-17-5 ethanol

Oral	LD50	7,060 mg/kg (rat)
------	------	-------------------

Inhalative	LC50/4 h	20,000 mg/l (rat)
------------	----------	-------------------

- **Primary irritant effect:**
- **On the skin:** Based on available data, the classification criteria are not met.
- **On the eyes:** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitization** Based on available data, the classification criteria are not met.
- **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

· **12.1 Toxicity**

· **Aquatic toxicity:**

CAS: 64-17-5 ethanol

EC50 72h	275 mg/l (Algae)
----------	------------------

EC50 48h	12,340 mg/l (daphnia)
----------	-----------------------

EC50 96h	13,000 mg/l (fish)
----------	--------------------

- **12.2 Persistence and degradability :** No further relevant information available.
- **12.3 Bioaccumulative potential :** No further relevant information available.
- **12.4 Mobility in soil :** No further relevant information available.

(Contd. on page 6)



Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

Trade name: DETARTRINE 150 Z

(Contd. of page 5)

- **Additional ecological information:**
- **General notes:**
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water.
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **12.5 Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **12.6 Other adverse effects :** No further relevant information available.

SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation :** Smaller quantities can be disposed of with household waste.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

SECTION 14: Transport information

- | | |
|---|---|
| · 14.1 UN-Number
· ADR, ADN, IMDG, IATA | Not classified as hazardous for transport |
| · 14.2 UN proper shipping name
· ADR, ADN, IMDG, IATA | Not classified as hazardous for transport |
| · 14.3 Transport hazard class(es)
· ADR, ADN, IMDG, IATA
· Class | Not classified as hazardous for transport |
| · 14.4 Packing group
· ADR, IMDG, IATA | Not classified as hazardous for transport |
| · 14.5 Environmental hazards:
· Marine pollutant: | No |
| · 14.6 Special precautions for user | Not applicable. |
| · 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code | Not applicable. |
| · UN "Model Regulation": | Not classified as hazardous for transport |

SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

GB

(Contd. on page 7)

Safety data sheet
According to 1907/2006/EC, Article 31

Printing date 02.12.2017

Version number 3

Revision: 02.12.2017

Trade name: DETARTRINE 150 Z

(Contd. of page 6)

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Relevant phrases

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H319 Causes serious eye irritation.

H332 Harmful if inhaled.

H335 May cause respiratory irritation.

· Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

VOC: Volatile Organic Compounds (USA, EU)

DNEL: Derived No-Effect Level (REACH)

PNEC: Predicted No-Effect Concentration (REACH)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

SVHC: Substances of Very High Concern

vPvB: very Persistent and very Bioaccumulative

Flam. Liq. 2: Flammable liquids – Category 2

Acute Tox. 4: Acute toxicity – Category 4

Skin Irrit. 2: Skin corrosion/irritation – Category 2

Eye Irrit. 2: Serious eye damage/eye irritation – Category 2

STOT SE 3: Specific target organ toxicity (single exposure) – Category 3