



PLASTALGIN ALGINATE

Quick - Easy - Clean.

Makes a lasting impression.

Plastalgin is a multi-purpose, dust-free alginate impression material that works quickly and easily. Plastalgin creates a mixed creamy homogenous consistency without excessive flow.

Plastalgin is ideal for:

- Case study materials
- Orthodontic models
- Opposing models
- Mouthguards
- Bleaching trays
- Provisionals

AVAILABLE IN 1 LB. BAGS:

Plastalgin Quick Set -
Fruit Flavor

Plastalgin Normal Set -
Fruit Flavor



Making Impressions Quick and Easy...

Plastalgin is a highly versatile alginate impression material that is simple to use, easy to handle, and available in 2 formulations. Whether you choose Plastalgin Normal Set (fruit flavored) or Plastalgin Quick Set (fruit flavored), both products are naturally hydrophilic and deliver precise impressions with an excellent reproduction of marginal details.

Making impressions safe and dust-free...

Plastalgin's dust-free formulation is specifically designed to reduce the amount of airborne particles dispersed to exposed individuals, contributing to a healthier for you and your patients.

	Plastalgin Normal Set	Plastalgin Quick Set
Mixing Time:	45 seconds	30 seconds
Total Working Time:	2:15	1:15
Setting Time:	3:10	1:45
Viscosity:	38mm	36mm
Reproduction of Details:	20m	50m
Deformation When Compressed:	12.5%	12.5%

Ordering Information

ITEM NUMBER	PLASTALGIN TYPE
01-C1000	PLASTALGIN QUICK SET, 17OZ. - FRUIT FLAVORED
01-C1100	PLASTALGIN NORMAL SET, 17OZ. - FRUIT FLAVORED

Visit www.septodont.com for more information about Plastalgin Alginate and other Septodont products
800-872-8305 | septodont.com

